

GP DETERGENT

CONCENTRATED LIQUID DETERGENT

- Mild - ideal for hand use
- Leaves washed items sparkling and grease free
- Long lasting foam

USES

A detergent concentrate containing additives for efficient grease removal.

A thick, green liquid consisting of detergents, grease-dissolving surfactants, foam boosters, preservative and dye. GP detergent is pH-balanced for mildness to hands.

Designed for cleaning glassware, cutlery, crockery, utensils and worktops, GP Detergent is ideal for use on aluminium, glass, stainless steel, vinyl, plastic, ceramic, wood and enamel.


DIRECTIONS

DISHWASHING

Dilute 1 part GP Detergent with between 750 to 1000 parts warm water in a washing-up bowl. Add items to be washed and allow a few minutes to soak. Clean items using a sponge, brush or pad.

GENERAL CLEANING



Dilute 1 part GP Detergent with 300 parts cold or warm water.



DILUTIONS

| GP DETERGENT : WATER | |
|-------------------------------|---------------------------|
| Dishwashing 1:750 - 1:1000 | General Cleaning 1:300 |

Code: 418



PACK SIZES



5L